



Inspection Report

UNIVERSITY OF CALIFORNIA, IRVINE

Customer ID: **9193**

Certificate: **93-R-0434**

Site: 001

160 ALDRICH HALL

UNIVERSITY OF CALIFORNIA, IRVINE

Type: ROUTINE INSPECTION

Date: Jun-17-2009

IRVINE, CA 92697

No non-compliant items identified this inspection.

Prepared By:

MANUEL A ADVIENTO, D.V.M. USDA, APHIS, Animal Care

Date:

Title: VETERINARY MEDICAL OFFICER Inspector 5001

Jun-17-2009

Received By:

(b)(6), (b)(7)(c)

Date:

Title:

Jun-17-2009



Inspection Report

UNIVERSITY OF CALIFORNIA, IRVINE

Customer ID: **9193**

160 ALDRICH HALL

Certificate: **93-R-0434**

Site: 001

UNIVERSITY OF CALIFORNIA, IRVINE

IRVINE, CA 92697

Type: ROUTINE INSPECTION

Date: Nov-05-2008

No non-compliant items identified this inspection.

This is a transcribed copy of the original inspection report.

Prepared By:

MANUEL A ADVIENTO, V.M.O. USDA, APHIS, Animal Care

Date:

Title: VETERINARY MEDICAL OFFICER Inspector 5001

Mar-23-2009

Received By:

(b)(6), (b)(7)(c)

Date:

Title:

Mar-23-2009



Inspection Report

UNIVERSITY OF CALIFORNIA, IRVINE

Customer ID: **9193**

160 ALDRICH HALL

Certificate: **93-R-0434**

Site: 001

UNIVERSITY OF CALIFORNIA, IRVINE

IRVINE, CA 92697

Type: ROUTINE INSPECTION

Date: Mar-11-2008

No non-compliant items identified rhis inspection.

Prepared By:

MANUEL A ADVIENTO, V.M.O. USDA, APHIS, Animal Care

Date:

Title: VETERINARY MEDICAL OFFICER Inspector 5001

Mar-11-2008

Received By:

(b)(6), (b)(7)(c)

Date:

Title:

Mar-11-2008



Inspection Report

UNIVERSITY OF CALIFORNIA, IRVINE

Customer ID: **9193**

Certificate: **93-R-0434**

Site: 001

UNIVERSITY OF CALIFORNIA, IRVINE

160 ALDRICH HALL

Type: ROUTINE INSPECTION

IRVINE, CA 92697

Date: Jul-17-2007

3.56 (a) (3)

SANITATION.

(3) Primary enclosures for rabbits shall be sanitized by washing them with hot water (180 °F.) and soap or detergent as in a mechanical cage washer

-There is no record to indicate that 180 degrees is attained in their mechanical washer. Pathogenic microorganisms may survive in temperatures below 180 degrees. Temperature monitoring is in place but nothing was done starting July 1, 2007. Regular monitoring shall be done to ensure that hot water temperature of 180 degrees is attained in their mechanical washer. CORRECT BY AUGUST 17, 2007.

Prepared By:

MANUEL A ADVIENTO, D.V.M. USDA, APHIS, Animal Care

Date:

Title: VETERINARY MEDICAL OFFICER Inspector 5001

Jul-17-2007

Received By:

(b)(6), (b)(7)(c)

Date:

Title:

Jul-17-2007